

Cannot export more than 100 records at a time

Frequently Asked Question:

Cannot export more than 100 records at a time

Answer:

You can experience problems with the export from WIAT-II Scoring Assistant when a large number of records are selected at once. This has been seen with as few as 100 records. Follow the work around below to get around this.

- Export 50 or so records only.
- Exit PsychCorpCenter.
- Open PsychCorpCenter again.
- Export another 50, and repeat step 2.

Article ID: 183

Last updated: 13 Apr, 2010

Software Support -> — W — -> WIAT-II Scoring Assistant -> Cannot export more than 100 records at a time

<http://pearsonassessmentsupport.com/support/index.php?View=entry&EntryID=183>