

Error 106 during the move data process.

Error Description:

Error 106 during the move data process.

Resolution:

The CMS Scoring Assistant software installs to the local hard drive. Please ensure that you have administrator access to the hard drive for installing software. If the problem persists, contact the Technical Support Center by email at esupport@pearson.com for assistance.

Article ID: 1895

Last updated: 14 Apr, 2010

Software Support -> — C — -> Childrens Memory Scale -> Error 106 during the move data process.

<http://pearsonassessmentsupport.com/support/index.php?View=entry&EntryID=1895>