Install PC or PDA software separately

Frequently Asked Question:

Install PC or PDA software separately

Answer:

The following are command line switches that can be used after the normal setup.exe command. This will allow you to install either the desktop software only, or the PDA software only.

- Installs PDA software only: -CompleteInstallWithPDA (ex. D:\Setup.exe -CompleteInstallWithPDA)
- Installs Desktop software only: -CompleteInstallNoPDA (ex. D:\Setup.exe -CompleteInstallNoPDA)

Note: A space must separate Setup.exe and the switch.

The Sensory Profile is no longer available for purchase. Please update your software to the enhanced Sensory Profile Select Scoring Assistant.

Article ID: 1040

Last updated: 01 Jul, 2010

Software Support -> -S --> Sensory Profile Scoring Assistant -> Install PC or PDA software separately

https://pearsonassessmentsupport.com/support/index.php?View=entry&EntryID=1040