

# The system cannot open the device or file specified

## Error Description:

The system cannot open the device or file specified

## Resolution:

Please go to <http://www.pearsonassessmentsupport.com/support/index.php> and select your windows platform and then other products. Locate the NEPSY-II Scoring Assistant folder, and update both the PsychCorpCenter 1 3.3.3, and the Nepys-II SA 2.0.1 programs.

---

Article ID: 1463

Last updated: 02 Jul, 2010

Software Support -> — N — -> NEPSY-II Scoring Assistant -> The system cannot open the device or file specified

<https://pearsonassessmentsupport.com/support/index.php?View=entry&EntryID=1463>