

## WPPSI-III ReadMe 2.0

### WPPSI-III Scoring Assistant Version 2.0

This file includes helpful information about the WPPSI-III Scoring Assistant, version 2.0. You may also view this file from the CD-ROM or through your Windows Explorer.

- 1 This application is compatible with windows Vista|2000|XP.
- 2 If you receive the following error message upon login, "TPCDTPicker Runtime error 2147208421 8004331b invalid window," applying this update will resolve the issue.
- 3 The preferred system settings for this software are 16-bit color, minimum of 800 x 600 display area, and small fonts. To check your settings, select Start | Settings | Control Panel. Click the 'Display' icon and select the 'Settings' tab. Click 'OK' to save any changes.
- 4 If you have version 1.0 of the WISC-III--WIAT-II Scoring Assistant on your desktop, you must install the WISC-III--WIAT-II Scoring Assistant update, version 1.1.2 after installing the WPPSI-III Scoring Assistant, but prior to running the application. You may download this update by going to <http://www.psychcorp.com> and click on the Software Support link.

For additional tips and information, read the online help after you have installed the WPPSI-III Scoring Assistant.