

Error 53 location PsychCorpCenter process new products tpcasmbengine.cls line 220

Error Description:

Error 53 location PsychCorpCenter process new products tpcasmbengine.cls line 220

Resolution:

This error occurs when you install the WISC IV SA version 1.0, 1.0.1 and 1.0.3 after updating your WAIS-III/WMS-III/WIAT II Scoring Assistant to version 1.0.5. To correct, please follow the following steps.

1. Use My Computer or Windows Explorer to locate C:\Program Files\PsychCorpCenter\Data.
2. Copy the file as _data.mdb on the desktop and rename it as _data.bkp.
3. Uninstall all previously installed PsychCorpCenter products.
4. Delete the folder PsychCorpCenter.
5. Empty the Recycle bin.
6. Install all PsychCorpCenter products back on to the machine, taking mind to the install order.

Open the software and do a restore of the backed up database on the desktop.

As another option to correct the issue, the newest PsychCorpCenter product and platform updates can be applied from our website www.pearsonassessmentsupport.com

Article ID: 94

Last updated: 30 Jun, 2010

Software Support -> — W — -> WISC-IV Scoring Assistant -> Error 53 location PsychCorpCenter process new products tpcasmbengine.cls line 220

<https://pearsonassessmentsupport.com/support/index.php?View=entry&EntryID=94>