

Run-time error 339 Component SFRMA35.OCX or one of its dependencies not correctly registered

Error Description:

Run-time error 339 Component SFRMA35.OCX or one of its dependencies not correctly registered

Runtime error - a file is missing or invalid

SPDTP Error -2147467259 (80004005) could not delete from specified folder

Runtime 2147467259(80004005) at logon

SP DDT: Caution Runtime error 94, invalid use of null

Resolution:

You can uninstall the software, and then reinstall the software without installing the PDA component unless you plan to use the PDA module. Please close background applications before installing the software.

Windows NT/2000/XP users must have the ability to read and write to the program folder. Some users may require additional rights.

The Sensory Profile is no longer available for purchase. Please update your software to the enhanced Sensory Profile Select Scoring Assistant.

Article ID: 994

Last updated: 01 Jul, 2010

Software Support -> — S — -> Sensory Profile Scoring Assistant -> Run-time error 339 Component SFRMA35.OCX or one of its dependencies not correctly registered

<https://pearsonassessmentssupport.com/support/index.php?View=entry&EntryID=994>