

Test Components and Equipment

SCAN-C: 3 kit components

- Manual
- 25 Record Forms
- audio compact disc (CD)

You will need to provide (purchase yourself)

- a CD player with track display or a two-channel audiometer
- two sets of stereo headphones: one set for the child one set for you to monitor the auditory stimuli
- a stereo Y-adaptor if the CD player does not have two output jacks
- SCAN-C: 3 can also be administered using a computer with a headphone jack into which a stereo Y-adaptor can be inserted.

For Additional information please click on the link below

Article ID: 3608

Last updated: 01 Jul, 2010

Software Support -> — S — -> SCAN-3: C -> Test Components and Equipment

<https://pearsonassessmentsupport.com/support/index.php?View=entry&EntryID=3608>